

Product Datasheet

CD282 Recombinant Protein (orb1495133)

Description	CD282 Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, 4M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	ESSNQASLSCDRNGICKGSSGSLNSIPSGLTEAVKSLDLSNNRITYISNSDLQRCV NLQALVLTNSGINTIEEDSFSSLGSLHLDLSYNYLSNLSSSWFKPLSSLTFLNLLG NPYKTLGETSLFSLTKLQILRVGNMDTFTKIQRKDFAGLTFLEELEIDASDLQSYE PKSLKSIQNVSHLILHMKQHILLLEIFVDVTSSVECLELRDLDLDTFHFSELSTGET NSLIKKFTFRNVKITDESLFQVMKLLNQISGLLELEFDDCTLNGVGNFRASDNR VIDPGKVETLTIRRLHIPRFYLFYDLSTLYSLTERVKRITVENSKVFLVPCLLSQHLK SLEYLDLSENLMVEEYLKNSACEDAWPSLQTLILRQNHLSLEKTGETLLTLKNLT NIDISKNSFHSMPEQCWPEKMKYLNLSSTRIHSVTGCIPKLEILDVSNNNLNLF SLNLPQLKELYISRNKMLTLPDASLLPMLLVLKISRNAITTFKSKEQLDSFHLLKTLLEA GGNNFICSCFEFLSFTQEQALAKVLIDWPANYLDCSPSHVVRGQVQDVRLSVSE CHRT
MW	~68kDa
Uniprot ID	O60603
Solubility (25°C)	PBS, 4M Urea, PH7.4
Expiration Date	6 months from date of receipt.