

Product Datasheet

CD62E Recombinant Protein (orb1495147)

Description	CD62E Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, 4M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	MWSYNTSTEAMTYDEASAYCQQRythlVAIQNKEEIEYLNsilSYSPSYWIGIRK VNNVWVWVGTQKPLTEEAKNWAPGEPNNRQKDEDCVEIYIKREKDVGMWNDE RCSKKLALCYTAACNTSCSGHGECVETINNYTCKCDPGFSGLKCEQIVNCTAL ESPEHGSLVCSHPLGNFSYNSSCSISCDRGLPSSMETMQCMSSGEWSAPIAC NVVECDAVTNPANGFVECFQNPGSFPWNTTCTFDCEEGFELMGAQSLQCTSSG NWDNEKPTCKAVTCRAVRQPQNGSVRCSHSPAGEFTFKSSCNFTCEEGFMLQG PAQVECTTQGQWTQQIPVCEAFQCTALSNPERGYMNCLPSASGSFRYGSSCEFS CEQGFVLKGSKRLQCGPTGEWDNEKPTCEAVRCDAVHQPPKGLVRCAHSPIGE FTYKSSCAFSCEEGFELHGSTQLECTSQGQWTEEVPCQVVKSSSLAVPGKINM SCSGEPVFGTVCKFACPEGWTLNGSAARTCGATGHWSGLLPTCEAPTESNIP
MW	~62kDa
Uniprot ID	P16581
Solubility (25°C)	PBS, 4M Urea, PH7.4
Expiration Date	6 months from date of receipt.