

Product Datasheet

CD299 Recombinant Protein (orb1495152)

| | |
|--------------------------|---|
| Description | CD299 Recombinant Protein |
| Species/Host | E. coli |
| Preservatives | PBS, 4M Urea, pH 7.4 |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Tag | His-tag |
| Note | For research use only |
| Protein Sequence | <p>MQVSKVPSSLSQEQSEQDAIYQNLTLKAAVGELSEKSKLQEIYQELTQLKAAVG ELPEKSKLQEIYQELTRLKAAVGELPEKSKLQEIYQELTRLKAAVGELPEKSKLQEI YQELTRLKAAVGELPEKSKLQEIYQELTELKAAVGELPEKSKLQEIYQELTQLKAAV GELPDQSKQQIYQELTDLKTAFERLCRHCPKDWTFQGNCFMSNSQRNWH DSVTACQEVRAQLVVIKTAEEQNFLQLQTSRSNRFSWMGLSDLNQEGTWQWV DGSPLSPSFQRYWNSGEPNNSGNEDCAEFSGSGWVNDNRCDVDNYWICKPAA CFRDE</p> |
| MW | ~38kDa |
| Uniprot ID | Q9H2X3 |
| Solubility (25°C) | PBS, 4M Urea, PH7.4 |
| Expiration Date | 6 months from date of receipt. |