

Product Datasheet

CD23A Recombinant Protein (orb1495157)

Description	CD23A Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, 4M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	MDTTQSLKQLEERAARNVSQVSKNLESHHGDQMAQKSQSTQISQELEELRAEQ QRLKSQDLELSWNLNGLQADLSSFKSQELNERNEASDLLERLREEVTKLRMELQ VSSGFVCNTCPEKWINFQRKCYFYGKGTKQVHARYACDDMEGQLVSIHSPEE QDFLTKHASHTGSWIGLRNLDLKGFEIWDGSHVDYSNWAPGEPTSRSQGEDC VMMRGSGRWNDADFCDRKLGAWVCDRLATCTPPASEGSAESMGPDSPDPDG RLPTPSAPLHS
MW	~35kDa
Uniprot ID	P06734
Solubility (25°C)	PBS, 4M Urea, PH7.4
Expiration Date	6 months from date of receipt.