

Product Datasheet

CD218b Recombinant Protein (orb1495166)

Description	CD218b Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, 4M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	MFNISGCSTKLLWTYSTRSEEEFVLFCDLPEPQKSHFCHRNLSPKQVPEHLPF MGSNDLSDVQWYQQPSNGDPLEDIRKSYPHIIQDKCTLHFLTPGVNNSGSYICR PKMIKSPYDVACCVKMILEVKPQTNASCEYSASHKQDLLLGSTGSISCPSLSCQS DAQSPAVTWYKNGKLLSVERSNRIVVDEVYDYHQGTIVCDYTQSDTVSSWTVR AVVQVRTIVGDTKLKPDILDPEVDTLEVELGKPLTISCKARFGFERVFNPKWYIK DSDLEWEVSVPEAKSIKSTLKDDEIERNIILEKVTQRDLRRKFVCFVQNSIGNTTQS VQLKEKR
MW	~43kDa
Uniprot ID	O95256
Solubility (25°C)	PBS, 4M Urea, PH7.4
Expiration Date	6 months from date of receipt.