

Product Datasheet

GRF-pan polyclonal antibody (orb1495237)

Description	GRF-pan polyclonal antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of arabidopsis thaliana GRF-pan. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 31 kDa
Uniprot ID	P48348 , Q96299 , P42644 , Q96300 , P42643 , P48349 , P46077 , Q01525 , P42645
Dilution Range	WB (1/500 - 1/2000)
Expiration Date	12 months from date of receipt.