

Product Datasheet

COVID-19 Spike glycoprotein polyclonal antibody (orb1495253)

Description	COVID-19 Spike glycoprotein polyclonal antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of COVID-19 Spike glycoprotein. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 40 kDa
Uniprot ID	P0DTC2
Dilution Range	WB (1/1000 - 1/3000), ELISA: Use at an assay dependent dilution
Expiration Date	12 months from date of receipt.