



## Product Datasheet

# p70 S6 Kinase (Phospho-Ser411) Antibody (orb14953)


|                            |   |
|----------------------------|---|
| <b>Description</b>         | p70 S6 Kinase (Phospho-Ser411) Antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IF, IHC, WB   |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Recombinant   |
| <b>Uniprot ID</b>          | <b>P23443</b>   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Confocal immunofluorescent analysis of M...



Immunohistochemical analysis of formalin...



Western blot analysis of p70 S6 Kinase(P...