

Product Datasheet

Claudin 7 monoclonal antibody (orb1495357)

Description	Claudin 7 monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Claudin 7. The exact sequence is proprietary.
Preservatives	Mouse IgG2a. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 22 kDa
Uniprot ID	O95471
Dilution Range	WB (1/500 - 1/1000), IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.