

## Product Datasheet

### Claudin 7 monoclonal antibody (orb1495376)

<b>Description</b>	Claudin 7 monoclonal antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence within human Claudin 7. The exact sequence is proprietary.
<b>Preservatives</b>	Mouse IgG2a. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	~ 22 kDa
<b>Uniprot ID</b>	<a href="#">O95471</a>
<b>Dilution Range</b>	WB (1/500 - 1/1000), IHC (1/100 - 1/300)
<b>Expiration Date</b>	12 months from date of receipt.