

Product Datasheet

GST-tag polyclonal antibody (orb1495631)

Description	GST-tag polyclonal antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant full length protein.
Preservatives	Liquid in PBS, pH 7.3, 0.2% rAlbumin, and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.