

Product Datasheet

Myc-tag polyclonal antibody (orb1495632)

Description	Myc-tag polyclonal antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	KLH-conjugated synthetic peptide EQKLISEEDL.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB (1/1000 - 1/2000), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.