

## Product Datasheet

### HA-tag polyclonal antibody (orb1495634)

<b>Description</b>	HA-tag polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide YPYDVPDYA.
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB (1/1000 - 1/2000), IP (1/10 - 1/100)
<b>Expiration Date</b>	12 months from date of receipt.