


Product Datasheet


p62Dok (phospho-Tyr362) Antibody (orb14957)

Description	p62Dok (phospho-Tyr362) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q99704
Expiration Date	12 months from date of receipt.



Western blot analysis of p62Dok(phospho-

...



Immunohistochemical analysis of formalin...