

Product Datasheet

PKC gamma (phospho-T655) polyclonal antibody (orb1495826)

Description	PKC gamma (phospho-T655) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human PKC gamma with a phosphorylation site at T655. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 78 kDa
Uniprot ID	P05129
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.