

Product Datasheet

PRKD1/2/3 (phospho-S738/S742) polyclonal antibody (orb1495947)

Description	PRKD1/2/3 (phospho-S738/S742) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human PRKD1/2/3 (pS738/S742). The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 130 kDa
Uniprot ID	Q15139 , O94806 , Q9BZL6
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.