

Product Datasheet

MST1/2 (phospho-T183) polyclonal antibody (orb1495954)

| | |
|----------------------------|---|
| Description | MST1/2 (phospho-T183) polyclonal antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human MST1/2 (pT183). The exact sequence is proprietary. |
| Preservatives | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 60 kDa |
| Uniprot ID | Q13188 , Q13043 |
| Dilution Range | WB (1/500 - 1/1000), IH (1/100 - 1/200) |
| Expiration Date | 12 months from date of receipt. |