



Product Datasheet

eIF2 α (Phospho-Ser51) Antibody (orb14960)


Description	eIF2 α (Phospho-Ser51) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05198
Expiration Date	12 months from date of receipt.



Western blot analysis of eIF2 α (Phospho-S...



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of H...