

Product Datasheet

NF-κB p65(Phospho-Ser276) polyclonal antibody (orb1496029)

Description	NF-κB p65(Phospho-Ser276) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human NF-κB p65 around the phosphorylation site of Serine 276.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 60 kDa
Uniprot ID	Q04206
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.