

## Product Datasheet

### GSK3 $\beta$ (Ab-9) polyclonal antibody [KU] (orb1496046)

<b>Description</b>	GSK3 $\beta$ (Ab-9) polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Immunogen</b>	Synthetic peptide, corresponding to Human GSK3 $\beta$ .
<b>Preservatives</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 46 kDa
<b>Uniprot ID</b>	<a href="#">P49841</a>
<b>Dilution Range</b>	WB: 1:1000~1:2000
<b>Expiration Date</b>	12 months from date of receipt.