

Product Datasheet

A35R Sheep polyclonal antibody (orb1496103)

Description	A35R Sheep polyclonal antibody
Species/Host	Sheep
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant protein of MPXV A35R.
Preservatives	PBS, pH 7.4
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 20 kDa
Uniprot ID	Q80KX2
Dilution Range	ELISA: 1:1000~1:30000
Expiration Date	12 months from date of receipt.