

## Product Datasheet

### A29L Sheep polyclonal antibody (orb1496104)

<b>Description</b>	A29L Sheep polyclonal antibody
<b>Species/Host</b>	Sheep
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant protein of MPXV A29L.
<b>Preservatives</b>	PBS, pH 7.4
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 12 kDa
<b>Uniprot ID</b>	<a href="#">Q9YN60</a>
<b>Dilution Range</b>	ELISA: 1:1000~1:30000
<b>Expiration Date</b>	12 months from date of receipt.