

Product Datasheet

Rabbit TROP-2 / TACSTD2 Protein (orb1496137)

Description	Rabbit TROP-2, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Pro 279 (...)
Reactivity	Rabbit
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TROP-2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.5 kDa. The protein migrates as 40-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Gln 35 - Pro 279
Purity	95%
Source	Rabbit TROP-2, His Tag (orb1496137) is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Pro 279 (Accession # G1T6B6-1).
MW	29.5 kDa
Uniprot ID	G1T6B6
Expiration Date	6 months from date of receipt.