

## Product Datasheet

### Mouse IgG3 Fc Protein (orb1496144)

<b>Description</b>	Mouse IgG3 Fc Protein, Tag Free is expressed from human 293 cells (HEK293). It contains AA Pro 98 - ...
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IgG3 Fc
<b>Preservatives</b>	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	Native
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries no "tag". The protein has a calculated MW of 26.1 kDa. The protein migrates as 30-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	Pro 98 - Pro 327
<b>Purity</b>	90%
<b>Source</b>	Mouse IgG3 Fc Protein, Tag Free (orb1496144) is expressed from human 293 cells (HEK293). It contains AA Pro 98 - Pro 327 (Accession # P03987-2).
<b>MW</b>	26.1 kDa
<b>Uniprot ID</b>	<a href="#">P03987</a>
<b>Expiration Date</b>	6 months from date of receipt.