

Product Datasheet

Rabbit HGF R / c-MET Protein (orb1496163)

Description	Rabbit HGF R Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Me...
Reactivity	Rabbit
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	HGF R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 103.3 kDa. The protein migrates as 35-43 kDa and 85-95 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001164511.1
Purity	90%
Source	Rabbit HGF R Protein, His Tag (orb1496163) is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Met 933 (Accession # Q09YN5-1).
MW	103.3 kDa
Uniprot ID	Q09YN5
NCBI	NP_001164511.1
Expiration Date	6 months from date of receipt.