

Product Datasheet

Mouse ALPG / ALPPL2 Protein (orb1496169)

Description	Mouse ALPG, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 19 - Ser 502 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ALPG
Preservatives	20mM Tris,150mM NaCl,pH 7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 54.5 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Val 19 - Ser 502
Purity	95%
Source	Mouse ALPG, His Tag (orb1496169) is expressed from human 293 cells (HEK293). It contains AA Val 19 - Ser 502 (Accession # P24823-1).
MW	54.5 kDa
Uniprot ID	P24823
Expiration Date	6 months from date of receipt.