
Product Datasheet

Cynomolgus Coagulation factor X / F10 Protein, His Tag (orb1496179)

Description

Cynomolgus Coagulation factor X /

Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Coagulation factor X
Preservatives	25 mM MES,150 mM NaCl,pH6.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.7 kDa and 30.3 kDa. The protein migrates as 19 kDa,48-53 kDa and 65-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Leu 24 - Lys 488
Purity	90%
Source	Cynomolgus Coagulation factor X, His Tag (orb1496179) is expressed from human 293 cells (HEK293). It contains AA Leu 24 - Lys 488 (Accession # A0A2K5VHK1-1).
MW	15.7 kDa and 30.3 kDa
Uniprot ID	A0A2K5VHK1
Expiration Date	6 months from date of receipt.