

Product Datasheet

Rat CD28 Protein(BLI verified) (orb1496194)

Description	Rat CD28, His Tag is expressed from human 293 cells (HEK293). It contains AA Asn 20 - Leu 150 (Acces...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD28
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.9 kDa. The protein migrates as 30-37 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_037253.1
Purity	90%
Source	Rat CD28, His Tag (orb1496194) is expressed from human 293 cells (HEK293). It contains AA Asn 20 - Leu 150 (Accession # P31042-1).
MW	16.9 kDa
Uniprot ID	P31042
NCBI	NP_037253.1
Expiration Date	6 months from date of receipt.