

Product Datasheet

Human Plectin / PLTN Protein (orb1496199)

Description	Human Plectin, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 4377 - Ala 468...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Plectin
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 34.5 kDa. The protein migrates as 43-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_958782.1
Purity	90%
Source	Human Plectin, His Tag (orb1496199) is expressed from human 293 cells (HEK293). It contains AA Ala 4377 - Ala 4684 (Accession # Q15149-1).
MW	34.5 kDa
Uniprot ID	Q15149
NCBI	NP_958782.1
Expiration Date	6 months from date of receipt.