
Product Datasheet

Cynomolgus OSMR beta Protein (orb1496200)

Description

Cynomolgus OSMR, His Tag is expressed

contains AA Glu 28 - Ser 737...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

OSMR

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 82.8 kDa. The protein migrates as 105-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Cynomolgus OSMR, His Tag (orb1496200) is expressed from human 293 cells (HEK293). It contains AA Glu 28 - Ser 737 (Accession # A0A2K5UFW5-1).

MW

82.8 kDa

Uniprot ID

[A0A2K5UFW5](#)

Expiration Date

6 months from date of receipt.