

Product Datasheet

Cynomolgus Growth Hormone R (GHR) Protein (orb1496205)

Description

Cynomolgus Growth Hormone R, His Tag is

 Cynomolgus Growth Hormone R, His Tag EISA
 1.1 µg of Cynomolgus Growth Hormone R, His Tag per well


	contains AA Phe ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Growth Hormone R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 30.2 kDa. The protein migrates as 40-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Cynomolgus Growth Hormone R, His Tag (orb1496205) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - Tyr 264 (Accession # EHH62357.1).
MW	30.2 kDa
Uniprot ID	G8F4Z3
NCBI	EHH62357.1
Expiration Date	6 months from date of receipt.

 Immobilized
 Human
 Growth
 Hormone
 Protein...