

Product Datasheet

Cynomolgus CD2 / SRBC Protein (orb1496216)

Description	Cynomolgus CD2, His Tag is expressed from human 293 cells (HEK293). It contains AA Lys 25 - Asp 209 ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.3 kDa. The protein migrates as 32-36 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Lys 25 - Asp 209
Purity	95%
Source	Cynomolgus CD2, His Tag (orb1496216) is expressed from human 293 cells (HEK293). It contains AA Lys 25 - Asp 209 (Accession # Q6SZ61-1).
MW	23.3 kDa
Uniprot ID	Q6SZ61
Expiration Date	6 months from date of receipt.