

Product Datasheet

Human TRAC Protein (orb1496223)

Description	Human TRAC, His Tag is expressed from human 293 cells (HEK293). It contains AA Ile 1 - Ser 115 (Acce...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TRAC
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 14.8 kDa. The protein migrates as 25-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Ile 1 - Ser 115
Purity	90%
Source	Human TRAC, His Tag (orb1496223) is expressed from human 293 cells (HEK293). It contains AA Ile 1 - Ser 115 (Accession # P01848-1).
MW	14.8 kDa
Uniprot ID	P01848
Expiration Date	6 months from date of receipt.