
Product Datasheet

Biotinylated Human ALPG / ALPPL2 Protein (orb1496228)

Description

Biotinylated Human ALPG Protein,

cells (HEK293). It contains AA Ile...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

ALPG

Preservatives

50 mM Tris,150 mM NaCl,pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 56.2 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Biotinylated Human ALPG Protein, His, (orb1496228) is expressed from human 293 cells (HEK293). It contains AA Ile 20 - Asp 503 (Accession # P10696-1).

MW

56.2 kDa

Uniprot ID**P10696****Expiration Date**

6 months from date of receipt.