

Product Datasheet

Mouse IgG2a Kappa Isotype Control (orb1496229)

Description	Mouse IgG2a Kappa Isotype Control
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG control
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower. Please avoid repeated freeze-thaw cycles
Note	For research use only
Application notes	This antibody is suitable for use as a non-targeting isotype control in various in vitro and in vivo studies. It can also be used as a negative control in various applications such as ELISA, Western blot, immunofluorescence, immunohistochemistry, immunoprecipitation, and flow cytometry. Each laboratory should determine an optimum working titer for use in its particular application.
Isotype	Mouse IgG2a Mouse kappa
Clonality	Recombinant
Purity	95%
Expiration Date	12 months from date of receipt.