

Product Datasheet

Cynomolgus IL-10 Protein, His Tag (orb1496244)

Description Cynomolgus IL-10 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 19 ...

Reactivity Monkey

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target IL-10

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag N-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 20.6 kDa. The protein migrates as 21-22 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.


Purity 90%

Source Cynomolgus IL-10 Protein, His Tag (orb1496244) is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Asn 178 (Accession # P79338).

MW 20.6 kDa

Uniprot ID **P79338**

Expiration Date 6 months from date of receipt.



Immobilized
Cynomolgus
IL-10 R
alpha Pro...