

## Product Datasheet

### Human CEACAM-6 / CD66c Protein (orb1496249)

<b>Description</b>	Human CEACAM-6, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Lys 35 - Gly 320 (...)
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CEACAM-6
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 57.7 kDa. The protein migrates as 75-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	Lys 35 - Gly 320
<b>Purity</b>	95%
<b>Source</b>	Human CEACAM-6, Fc Tag (orb1496249) is expressed from human 293 cells (HEK293). It contains AA Lys 35 - Gly 320 (Accession # P40199-1).
<b>MW</b>	57.7 kDa
<b>Uniprot ID</b>	<b>P40199</b>
<b>Expiration Date</b>	6 months from date of receipt.