

---

## Product Datasheet

# Monkeypox virus (strain Zaire-96-I-16) L1R Protein (orb1496250)

**Description**

Monkeypox virus (strain Zaire-96-I-

human 293 cells (HEK293). It c...

**Reactivity**

Virus

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

L1R

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

N-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 19.7 kDa. The protein migrates as 18-19 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

NP\_536512.1

**Purity**

95%

**Source**

Monkeypox virus (strain Zaire-96-I-16) L1R, His Tag (orb1496250) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Gln 152 (Accession # Q8V4Z7).

**MW**

19.7 kDa

**Uniprot ID**[Q8V4Z7](#)**NCBI**[NP\\_536512.1](#)**Expiration Date**

6 months from date of receipt.