

Product Datasheet

Human IL-17 RB Protein (orb1496253)

Description

Human IL-17 RB, Fc Tag is expressed from

Arg 18 - Pro 292 (...)

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-17 RB

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH 7.5

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 56.7 kDa. The protein migrates as 65-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human IL-17 RB, Fc Tag (orb1496253) is expressed from human 293 cells (HEK293). It contains AA Arg 18 - Pro 292 (Accession # Q9NRM6-1 (EVI27) & Q9NRM6-1(IL17Rh1)).

MW

56.7 kDa

Uniprot ID
Q9NRM6
NCBI
NP_061195.2
Expiration Date

6 months from date of receipt.