

Product Datasheet

Human NT-4 / Neurotrophin-4 Protein (orb1496255)

Description	Human NT-4, Tag Free is expressed from E. coli cells. It contains AA Gly 81 - Ala 210 (Accession # P...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	NT-4
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 14.1 kDa. The protein migrates as 15 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	Gly 81 - Ala 210
Purity	95%
Source	Human NT-4, Tag Free (orb1496255) is expressed from E. coli cells. It contains AA Gly 81 - Ala 210 (Accession # P34130-1).
MW	14.1 kDa
Uniprot ID	P34130
Expiration Date	6 months from date of receipt.