

Product Datasheet

Human LYPD1 / PHTS protein (orb1496290)

Description	Human LYPD1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Gly 115 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LYPD1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 36.8 kDa . The protein migrates as 43 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Leu 21 - Gly 115
Purity	95%
Source	Human LYPD1, Fc Tag (orb1496290) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Gly 115 (Accession # Q8N2G4-1).
MW	36.8 kDa
Uniprot ID	Q8N2G4
Expiration Date	6 months from date of receipt.