

Product Datasheet

Human Alpha-Synuclein (A53T) Protein, His Tag (orb1496292)

Description	Human Alpha-Synuclein (A53T) Protein, His Tag
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Alpha-synuclein
Preservatives	100 mM NaAC, pH7.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.4 kDa. The protein migrates as 18-19 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	Met 1 - Ala 140
Purity	90%
Source	Human Alpha-Synuclein (A53T), His Tag (orb1496292) is expressed from E. coli cells. It contains AA Met 1 - Ala 140 (Accession # P37840-1 (A53T)).
MW	16.4 kDa
Uniprot ID	P37840
Expiration Date	6 months from date of receipt.