

Product Datasheet

Cynomolgus PD-1 / PDCD1 Protein (orb1496300)

Description	Cynomolgus PD-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 17.8 kDa. The protein migrates as 30-40 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Leu 25 - Gln 167
Purity	90%
Source	Cynomolgus PD-1, His Tag (orb1496300) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Accession # B0LAJ3).
MW	17.8 kDa
Uniprot ID	B0LAJ3
Expiration Date	6 months from date of receipt.