
Product Datasheet


Cynomolgus IL-2 R alpha / CD25 Protein (orb1496308)

Description

Cynomolgus IL-2 R alpha Protein, His Tag is

Cynomolgus IL-2 R alpha (CD25, His Tag) EA
 1.0 µg of Cynomolgus IL-2 R alpha (CD25, His Tag) EA

	contains AA ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-2 R alpha
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.3 kDa. The protein migrates as 52-57 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Cynomolgus IL-2 R alpha Protein, His Tag (orb1496308) is expressed from human 293 cells (HEK293). It contains AA Glu 22 - Glu 238 (Accession # H6WS54).
MW	26.3 kDa
Uniprot ID	H6WS54
Expiration Date	6 months from date of receipt.


 Immobilized
 Cynomolgus
 IL-2 R alpha
 / CD...