

Product Datasheet

Human TRDC Protein (orb1496315)

Description	Human TRDC, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 1 - Val 129 (Acce...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TRDC
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.4 kDa. The protein migrates as 27-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Ser 1 - Val 129
Purity	90%
Source	Human TRDC, His Tag (orb1496315) is expressed from human 293 cells (HEK293). It contains AA Ser 1 - Val 129 (Accession # B7Z8K6-1).
MW	16.4 kDa
Uniprot ID	B7Z8K6
Expiration Date	6 months from date of receipt.