

## Product Datasheet

### Human FXN ELISA Kit (orb1496419)

<b>Description</b>	Human FXN ELISA Kit
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Range</b>	12.5-800 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>Q16595</b>
<b>Sensitivity</b>	5.8 pg/mL
<b>Expiration Date</b>	Please enquire.

