

## Product Datasheet

### Zebrafish NE ELISA Kit (orb1496487)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Species/Host</b>	Zebrafish
<b>Reactivity</b>	Zebrafish
<b>Range</b>	78.13-5000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	26.1 pg/mL
<b>Expiration Date</b>	Please enquire.

