

## Product Datasheet

### Mouse Prok2 ELISA Kit (orb1496516)

<b>Description</b>	Mouse Prok2 ELISA Kit
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse
<b>Range</b>	0.32-20 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Uniprot ID</b>	<b>Q9QXU7</b>
<b>Sensitivity</b>	0.121 ng/mL
<b>Expiration Date</b>	Please enquire.

