

## Product Datasheet

### Mouse GPC3 ELISA Kit (orb1496525)

<b>Description</b>	Mouse GPC3 ELISA Kit
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse
<b>Range</b>	0.16-10 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates and other biological fluids
<b>Uniprot ID</b>	<a href="#">Q8CFZ4</a>
<b>Sensitivity</b>	0.058 ng/mL
<b>Expiration Date</b>	Please enquire.

