

## Product Datasheet

### Human HA ELISA Kit (orb1496533)

<b>Description</b>	Human HA ELISA Kit
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Range</b>	1.57-100 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	0.51 ng/mL
<b>Expiration Date</b>	Please enquire.

